



PRO ACTIV Reha-Technik GmbH
 Im Hofstätt 11
 D-72359 Dotternhausen
 Phone +49 (0)7427 9480-0
 Fax +49 (0)7427 9480-7025
 E-mail: info@proactiv-gmbh.de
 www.proactiv-gmbh.com



Offer date: _____

Order date: _____

according to GTC PRO ACTIV Reha-Technik GmbH

Consultant: _____

* Required fields, please fill in completely!

Order address / customer:

Cust. no.:* _____

Name: _____

First name: _____

Date of birth: _____

Street: _____

Postcode / city: _____

Phone:* (home) _____

Phone: (work) _____

E-Mail: _____

Collection: Yes No

Biometric data:

Body height: _____

Weight: _____

Inside leg: _____

FULL-E

€ 4,450.42

Integrated battery 630 Wh

- Aluminium hydroformed frame
- Tapered head tube 1 1/8" - 1.5"
- Seat post 31.6 mm
- Hydraulic disc brakes
- 148 mm BOOST standard
- Aluminium hollow chamber rim
- Shimano EP8 mid-motor
(250 W, 85 Nm, weight: 2.6 kg)

Ready-to-ride weight: starting at 25.5 kg

Permitted total weight: 120 kg

Frame size:

S (43 cm) A9515000100-043

M (47 cm) A9515000100-047

Chassis:

Tour

A9515100100

Forks: Manitou Machete 120mm

Dampers: Manitou MARA Inline 150 mm

Trail

A9515100110 **€ 244.71**

Forks: Manitou Mezzler Expert 160 mm

Dampers: Manitou MARA Inline 150 mm

Enduro

A9515100120 **€ 428.24**

Forks: Manitou Mezzler Pro 180 mm

Dampers: Manitou MARA Pro TR 180 mm

Wheel size

27.5" Wheels

A9500200000-584

**Combined wheel size:
Front 29" & rear 27,5"**

A9515200000 **€ 69.08**

Schwalbe Smart Sam 27.5x2.6"

A9500200050
(65-584)
Performance Line
RaceGuard 50



Schwalbe Rocket Ron 27.5x2.8"

A9500200051
(70-584)
Evolution Line
E25



Schwalbe Hans Dampf 27.5x2.8"

A9500200053
(70-584)
Evolution Line
E25



**Schwalbe Smart Sam 27.5x2.6"
Bronze Sidewall**

A9500200150 **€ 38.38**
(65-584)
Double Defence, RaceGuard E50
ADDIX



29" Wheels

A9500200010-622

Schwalbe Smart Sam Plus 29x2.25"

A9500200130
(57-622)
Performance Line
GreenGuard E50



Schwalbe Hans Dampf 29x2.35"

A9500200136
(65-622)
Evolution Line
E25



Dérailleur

Shimano Deore 1x10

A9515200100
Linkglide
Cassette: 11-43 teeth

Shimano XT 1x11

A9515200110 **€ 135.04**
Linkglide
Cassette: 11-50 teeth

Brake system

**Hydraulic disk brake system
Shimano Deore M6100
203 mm / 180 mm**

A9515300100

**Powerful hydraulic disk brake system
Shimano XT M8100 203 mm / 180 mm**

A9515300200 **€ 162.02**

**Disk brake system XT M 8120
Brake disk XT Ice Tech
203 mm / 180 mm**

A9515300300 **€ 241.18**

Saddle support

Aluminium saddle support 31.6 mm

A9500800300
Length: 350 mm



Lowerable saddle support 31.6 mm frame height 43

A9515800010 € 157.31
Remote on the handlebar
Length: 445 mm
Stroke: 150 mm



Lowerable saddle support 31.6 mm frame height 47

A9515800011 € 157.31
Remote on the handlebar
Length: 495 mm
Stroke: 175 mm



Saddle

Velo Fit Traveller saddle

A9500800010
Size M
Width: 15 cm; length: 27 cm



Selle Royal SCIENTIA saddle

A9500800*_____ € 66.39

Athletic (45°)

- Size 1 (small) (*103)
- Size 2 (medium) (*104)
- Size 3 (large) (*105)

Moderate (60°)

- Size 1 (small) (*106)
- Size 2 (medium) (*107)
- Size 3 (large) (*108)

Relaxed (90°)

- Size 1 (small) (*109)
- Size 2 (medium) (*110)
- Size 3 (large) (*111)



Front end

Handlebar stem 50 mm

Length: 50 mm
A9515400010-050
Clamping: Oversize 31.8 mm



Handlebar stem 70 mm

Length: 70 mm
A9515400010-070
Clamping: Oversize 31.8 mm



Ahead spacer 10 mm

A9500400500-011
2 pcs. included as standard
 A9500400500-010 € 0.84 per pc.
Additional installation: _____ pcs.



Handlebar

Sport Ergotec handlebar 720 mm

A9500400610
Width: 720 mm
Height: 15 mm
Grip bend: 9°
Oversize 31.8 mm



Sport Ergotec handlebar 780 mm

A9500400631
Width: 780 mm
Height: 15 mm
Grip bend: 9°
Oversize 31.8 mm



Grips

Sport grips black / grey

A9500400700
Shock-absorbing gel grip surface,
anti-twist clamping on the handlebar



Gel Ergo Wing grips

A9500400710 € 21.26
Ergonomic gel grip surface,
anti-twist clamping on the handlebar



Gel Ergo Wing Comfort grips

A9500400720 € 26.81
Large, ergonomically shaped wrist rest,
anti-twist clamping on the handlebar



Pedals

Offroad pedals

A9500500300
CR-MO axle,
with reflector (K-mark)
Platform with jagged profile



All-round Shimano PD-M324 pedals

A9500500350 € 50.67
CR-MO axle
One side with click system (SPD) and
adjustable release force,
other side with jagged profile



Pedals Shimano PD-GR500

A9500500340 € 63.40
CR-MO axle
Platform pedal with interchangeable screw pins
for optimum grip and control



Accessories

SIGMA LED lighting set

A9500900100 € 60.08
LED Lighting set compliant with road traffic
regulations
Consisting of front light and rear lights,
chargeable via Micro USB socket,
incl. charging cable.



SKS rear mudguard

A9500700052 € 25.13
detachable via quick-lock system




Stand with screwed mount

A9515900010 € 69.08
Colour: black

Bottle holder

A9500900750 € 11.26
incl. drinking bottle (600 ml)



Magnetic mechanical bottle holder system

A9500900751 € 39.33
incl. holder & drinking bottle (590 ml)



Bolted luggage carrier with mudguard

A9515700010 € 332.10
Only for rear wheel 27.5"
Loading up to 25 kg

IVERA Cable lock

A9500900720 € 36.72
Length: 85 cm
Security level: 6



IVERA Steel O-Flex cable lock








A9500900723 € 52.35
Length: 110 cm
Security level: 7



uGrip Bordo 5700/80 bicycle lock, black

A9500900700 € 51.51
Length: 80 cm
incl. transport bag
Security level: 7



<p>Bicycle lock Bordo 6000K 120 cm</p> <p><input type="checkbox"/> A9500900772 € 100.80</p> <p>Length: 120 cm incl. transport bag Security Level: 10</p> 	<p>Service kit</p> <p>Micro 200 Expandable saddle bag Incl. "SUPERSHORT" mini-pump, SUPER B TB-1118 tyre repair kit, M-WAVE Little multi-tool, SUPER B TB-5566 tyre lever hose 27.5" or 29"</p> <p><input type="checkbox"/> A9500900760-001 € 79.83 with hose 27.5"</p> <p><input type="checkbox"/> A9500900760-002 € 77.31 with hose 29"</p> 	<p>Micro 200 Expandable saddle bag</p> <p><input type="checkbox"/> A9500900120 € 33.53</p> <p>Dimensions: 16 cm x 19 cm x 12 cm Weight 175 g Volume: 1.5 l Load / cargo load: max. 1 kg Incl. saddle adapter</p> 
<p>"SUPERSHORT" mini-pump</p> <p><input type="checkbox"/> A9500900748 € 21.68</p> <p>Valve connection: AV, SV, DV with telescope function ergonomic T-grip Material: Aluminium / plastic Weight: 103 g</p> 	<p>SUPER B TB-1118 tyre repair kit</p> <p><input type="checkbox"/> A9500900746 € 2.18</p> <p>Set with 6 pcs. Adhesive tyre patches with sandpaper</p> 	<p>M-WAVE Little multi-tool</p> <p><input type="checkbox"/> A9500900747 € 18.15</p> <p>15 functions: 2/2.5/3/4/5/6/8 mm Allen, 1 flat head & 2 Phillips screwdrivers, Torx25, chain tool, 4 spoke spanners</p> 
<p>SUPER B TB-5566 tyre lever</p> <p><input type="checkbox"/> A9500900745 € 7.98</p> <p>Set with 3 pcs. Plastic with sturdy steel core with 2 hooks for hooking onto the spokes</p> 	<p>Hose 27.5"</p> <p><input type="checkbox"/> A9500900627 € 10.84 (54/75-584) AV 21F</p>	<p>Hose 29"</p> <p><input type="checkbox"/> A9500900266 € 8.32 (40/62-622/635) AV 19 AGV 40 mm</p>

Colour design – frame and attachment parts

For the following components, a colour can be chosen from the information sheet “Pedelec Colour Design”:
RAL powder paint or “Dormant” powder paint or “Liquid coating” custom paint

Frame colour:

RAL powder paint and “Dormant” powder paint:

Please enter the last 3 digits of the colour:

A9500100170 - in words (colour): _____

“Liquid coating” custom paint:

Please enter the last 2 digits of the custom colour:

€ 161.68

A9500100170 - **5** in words (colour): _____

Mudguard colour:

RAL powder paint or “Dormant” powder paint:

A9500700170

Please enter the last 3 digits of the powder paint:

A9500100170 - in words (colour): _____

“Liquid coating” custom paint:

A9500700172

Please enter the last 2 digits of the custom colour:

€ 111.85

A9500100170 - **5** in words (colour): _____

Notes and remarks:

Leasing

Yes

Still being discussed

No

Leasing companies: _____ (please fill in the “Pedelec leasing form”!)

I hereby confirm the information provided above and agree to receiving a written offer without obligation.

 Location, date End customer signature (offer)

 Location, date Consultant signature (offer)

 Location, date End customer signature (order)

 Location, date Consultant signature (order)

Info sheet: Colour Designs Pedelecs

Colour Designs - Powder Coatings

Information on current colour designs and sample images at: www.proactiv-gmbh.com/colours/

Powder Coatings

white RAL 9016 <input checked="" type="checkbox"/> gloss A9500100170-023	anthracite RAL 7016 <input checked="" type="checkbox"/> gloss A9500100170-001 <input checked="" type="checkbox"/> matt A9500100170-101	black RAL 9005 <input checked="" type="checkbox"/> gloss A9500100170-003 <input checked="" type="checkbox"/> matt A9500100170-019
pure orange RAL 2004 <input checked="" type="checkbox"/> gloss A9500100170-022 <input checked="" type="checkbox"/> matt A9500100170-122	yellow RAL 1023 <input checked="" type="checkbox"/> gloss A9500100170-046 <input checked="" type="checkbox"/> matt A9500100170-146	red RAL 3001 <input checked="" type="checkbox"/> gloss A9500100170-047 <input checked="" type="checkbox"/> matt A9500100170-147
„Dormant“ Powder Coatings		
sky blue <input checked="" type="checkbox"/> gloss A9500100170-017	electric green <input checked="" type="checkbox"/> gloss A9500100170-015	violet <input checked="" type="checkbox"/> gloss A9500100170-044
fire red <input checked="" type="checkbox"/> gloss A9500100170-018		
Powder Coatings “Trend colours”		
athletic grey <input checked="" type="checkbox"/> gloss A9500100170-303 € 30.59 <input checked="" type="checkbox"/> matt A9500100170-304 € 30.59	muddy grey <input checked="" type="checkbox"/> gloss A9500100170-305 € 30.59 <input checked="" type="checkbox"/> matt A9500100170-306 € 30.59	lagoon blue <input checked="" type="checkbox"/> gloss A9500100170-301 € 30.59 <input checked="" type="checkbox"/> matt A9500100170-302 € 30.59

Colour Designs - Custom Colours ‘Liquid Coating’

€ 161.68

storm (ice blue) <input type="checkbox"/> gloss A9500100170-500 <input type="checkbox"/> matt A9500100170-501	coffee <input type="checkbox"/> gloss A9500100170-520 <input type="checkbox"/> matt A9500100170-521	prairie (rose gold) <input type="checkbox"/> gloss A9500100170-504 <input type="checkbox"/> matt A9500100170-505
blue ocean <input type="checkbox"/> gloss A9500100170-522 <input type="checkbox"/> matt A9500100170-523	mint <input type="checkbox"/> gloss A9500100170-528 <input type="checkbox"/> matt A9500100170-529	