



EXPERIENCE MOBILITY WITHOUT COMPROMISES.

The **SPEEDY F2** is your individual rigid frame wheelchair with welded frame for unbeatable handling characteristics on flat or uneven paths. Each frame is unique according to your measurements and wishes – perfect for experienced wheelchair users. Select 14 chassis parameters such as knee angle, frame retraction, wheelbase and lower leg width at the top and bottom. The SPEEDY F2 is incredibly versatile with its individual frame designs, the different back versions and the various options.

SPEEDY F2



#FreeConception

Everyday companion and sports equipment in one. Whether it's dancing, paragliding, tennis, crossfit or world-class table tennis, the SPEEDY F2 can be flexibly adapted to your activities.



#Individuality

The dimensions are so flexible that the frame can be designed for people of short stature as well as for people over 2 m tall or leg amputees.



#Reliability

The SPEEDY F2 impresses with a sturdy and durable frame that ensures maximum mobility and is compatible with various additional drives and adaptations.



#Independence

Thanks to the easily foldable backrest and tool-free removal of assemblies, the SPEEDY F2 is easy to transport.



PERFECTION IN EVERY DETAIL.

Your **SPEEDY F2** is manufactured exactly according to your needs. You can choose from a wide range of options – from frame dimensions to footrests, brakes and push handles to back and seat systems, clothing guards, armrests and drive and caster wheels. Your wheelchair, your decision.

UNIQUE, SELECTABLE OPTIONS FROM PRO ACTIV:



The 7-way angle-adjustable and foldable backrest offers maximum flexibility for the ideal sitting position.



Transform your wheelchair into a piece of sports equipment in no time at all with the tool-free removable centre anti-tipping support.



You can choose between a permanently welded or removable wheelbase extension for more traction with a handbike or traction device.



With the Heavy Duty equipment, your SPEEDY F2 is designed for extreme loads and is suitable for user weights up to 200 kg.




TECHNICAL SPECIFICATIONS

Seat width:	25 - 55 cm (in 1 cm steps)
Seat depth:	23 - 55 cm (in 1 cm steps)
Back height:	20 - 60 cm (in 1 cm steps)
Lower leg length:	> 20 cm (in 1 cm steps)
Rear seat height:	35 - 50 cm (in 1 cm steps)
Front seat height:	36 - 58 cm (in 1 cm steps)
Tipping point:	-5 - 17 cm (in 1 cm steps)
Total product weight:	Starting at 8.5 kg
Max. standard load weight:	Up to 120 kg
Max. load weight for Heavy Duty:	Up to 200 kg appropriate individual dimensions upon request

Data subject to design changes.

FOLLOW US ON SOCIAL MEDIA



PROACTIV[®]

PRO ACTIV Reha-Technik GmbH

Im Hofstätt 11
72359 Dotternhausen
info@proactiv-gmbh.de

www.proactiv-gmbh.com

Rev. 09/2024

Subject to technical modifications, printing errors & colour deviations.