



BUDDY classic Active easy running wheelchair

Maximum user weight limit 100 kg

A9041100000 € 3,830.00

Quotation date: _____

Consultation was carried out by
Staff at the medical supplies shop: _____

Order date: _____

according to GTC PRO ACTIV Reha-Technik GmbH

Medical suppliers / company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Shipping address: same as company address

Other (see below)

Name: _____
Street: _____
Postcode / city: _____

Customer name / commission:


Surname, forename: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
D.O.B: _____
Paying organisation: _____

Invoicing address: same as company address

other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Body size:



Body height: _____ cm

Weight: _____ cm



Lower leg length: _____ cm

(anatomic dimension: hollow of the knee to the foot sole)

Indication: _____

This information is necessary for checking. If the information is missing, a check of the relevant items cannot be made!

Abducted frame shape, 2.5 cm/side



Lower leg length **USL**
(wheelchair dimension in consideration
of the seat cushion height)

Seat width **SB**

<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	<input type="checkbox"/> 31
<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38					
The following dimensions are optionally available subject to surcharge:				<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46

Size M (Frame length 34 cm (to the bend) and frame height 21 cm)

<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	<input type="checkbox"/> 31
-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------

Seat depth **ST**

Size L (Frame length 38 cm (to the bend) and frame height 21 cm)

<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35
-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------

Size XL (Frame length 43 cm (to the bend) and frame height 21 cm)

<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40
-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------

Back height **RH**

<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29
<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39
<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45				
Optional with size XL:			<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48				

Lower leg length **USL**

<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30
<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41
<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44								
Optional with size L + XL:			<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47					

Seat height back **SHh**

for 20" drive wheel						
<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37
for 22" drive wheel						
<input type="checkbox"/> 33,5	<input type="checkbox"/> 34,5	<input type="checkbox"/> 35,5	<input type="checkbox"/> 36,5	<input checked="" type="checkbox"/> 37,5	<input type="checkbox"/> 38,5	<input type="checkbox"/> 39,5
for 24" drive wheel						
<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input checked="" type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42

Seat height front **SHv**
(for 5" caster wheel)

<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47
-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------

Tipping point **KP**

(front side rear tube to
axis middle drive wheel)

<input checked="" type="checkbox"/> 5	<input type="checkbox"/> 6	<input type="checkbox"/> 7	<input type="checkbox"/> 8	<input type="checkbox"/> 9	<input type="checkbox"/> 10
---------------------------------------	----------------------------	----------------------------	----------------------------	----------------------------	-----------------------------

Please cross on the appropriate dimension (all dimensions are specified in cm).

Underlined values are standard settings.

The values SHh and SHv may deviate a little from the specified values depending on the wheel size, wheel camber as well as the front and rear frame angle!

Frame colour - one-colour design

Information on current colour designs and example images at: www.proactiv-gmbh.de/farben.php

Frame colour: Selection from colour palette

Frame colour will be specified with order

A9000100170-000

Raw, blasted, not coated

A9000100170-020

Note: In combination with moisture e.g. hand perspiration, can affect dark aluminum discoloration.

Powder colours

White RAL 9016

gloss A9000100170-023

matt A9000100170-123 € 42.10

Brilliant silver

gloss A9000100170-011

Anthracite RAL 7016

gloss A9000100170-001

matt A9000100170-101 € 42.10

Black RAL 9005

gloss A9000100170-003

matt A9000100170-019 € 42.10

Pure orange RAL 2004

gloss A9000100170-022

matt A9000100170-122 € 42.10

Pink magenta RAL 4010

gloss A9000100170-037

matt A9000100170-137 € 42.10

Pastel blue RAL 5024

gloss A9000100170-036

matt A9000100170-136 € 42.10

Lemon green

gloss A9000100170-005

matt A9000100170-105 € 42.10

Yellow RAL 1023

gloss A9000100170-046

matt A9000100170-146 € 42.10

Red RAL 3001

gloss A9000100170-047

matt A9000100170-147 € 42.10

Lagoon blue

gloss A9000100170-052

matt A9000100170-152 € 42.10

Athletic grey

gloss A9000100170-050

matt A9000100170-150 € 42.10

Muddy grey

gloss A9000100170-051

matt A9000100170-151 € 42.10

Powder colours „Dormant“

Deep black effect silver

gloss A9000100170-007

Sky blue

gloss A9000100170-017

Navy blue

gloss A9000100170-045

matt A9000100170-145 € 42.10

Electric green

gloss A9000100170-015

matt A9000100170-115 € 42.10

Violet

gloss A9000100170-044

Fire red

gloss A9000100170-018

matt A9000100170-118 € 42.10

Frame colour - One-colour design special colours “liquid coating”

€ 339.40

Storm (ice blue)

gloss A9000100170-500

matt A9000100170-501

Coffee

gloss A9000100170-520

matt A9000100170-521

Prairie (rose gold)

gloss A9000100170-504

matt A9000100170-505

Pearl (mother of pearl-white)

gloss A9000100170-524

matt A9000100170-525

Berry

gloss A9000100170-526

matt A9000100170-527

Blue ocean

gloss A9000100170-522

matt A9000100170-523

Emerald

gloss A9000100170-514

matt A9000100170-515

Mint

gloss A9000100170-528

matt A9000100170-529

Chili

gloss A9000100170-516

matt A9000100170-517

Frame colour - Two-colour design

Note: In the two-colour design, both colours are the same gloss level as the design colour.

Two-colour design with one base colour (on the rear frame) and one design colour (on the front frame) € 260.70

Please select a color from below:

Base colour (on the rear frame)

White RAL 9016

A9000100170-200

Brilliant silver

A9000100170-201

Anthracite RAL 7016

A9000100170-202

Black RAL 9005

A9000100170-203

Pastel blue RAL 5024

A9000100170-206

Lemon green

A9000100170-207

Yellow RAL 1023

A9000100170-208

Lagoon blue

A9000100170-212

Athletic grey

A9000100170-210

Muddy grey

A9000100170-211

Design colour (on the front frame)

Black RAL 9005

gloss A9000100170-259

matt A9000100170-260

Lemon green

gloss A9000100170-250

matt A9000100170-251

Chili

gloss A9000100170-253

matt A9000100170-254

Pure orange RAL 2004

gloss A9000100170-255

matt A9000100170-256

Coffee

gloss A9000100170-278

matt A9000100170-279

Mint

gloss A9000100170-286

matt A9000100170-287

Emerald

gloss A9000100170-274

matt A9000100170-275

Storm (ice blue)

gloss A9000100170-264

matt A9000100170-265

Prairie (rose gold)

gloss A9000100170-266

matt A9000100170-267

Pearl (mother of pearl-white)

gloss A9000100170-282

matt A9000100170-283

Berry

gloss A9000100170-284

matt A9000100170-285

Blue ocean

gloss A9000100170-280

matt A9000100170-281

Colour - Attachment parts

Note: If not specified, the caster fork and bearing block will be mounted in black.

Drive wheel: aluminum rim, hub and spokes

Black A9000100171-003

Caster fork and bearing block

Silver A9000100173-011

Black A9000100173-003

Colour design of attachment parts according to customer request:

Drive wheel: Aluminium rim and design inlet¹

Frame colour A9000100171-999 € 242.30

Select from one-colour design (different from the frame color) A9000100171-998 € 280.50

Colour: _____

A9000100170-

Caster wheel: 5" aluminium rim² (3-spoke design)

Frame colour A9000100174-999 € 73.00

Select from one-colour design (different from the frame color) A9000100174-998 € 84.60

Colour: _____

A9000100170-

Special colour requests (frame or attachment parts)

Frame colour for single purchasing

Selection from RAL colour palette

A9000100170-100 € 402.30

Gloss level

gloss

matt

colour name

RAL number

Note: attachment parts are **not** possible in design colors

¹ Only in conjunction with Liteline hub (with drive wheel aluminum rim with 18 or 36 spokes, single rim)

² Only in conjunction with caster wheel 5" aluminum rim 3-spoke design good shock absorption A900050121*, as well as with suspension fork 5-spoke design A9000500160*

Frames

Custom-made seat width dimension 39 cm - 46 cm


A9041100060 € 297.00

with frame reinforcement, double cross-tube bolted, side panel screw connection via the reinforcement inlet

Caster wheel bearing block

Base plate for bearing block welded

A9041500150



Base plate for bearing block clamped

A9041500600-K* € 114.00



Push handles

Safety push handles back-positioned

Length: 45 cm
A9041401075-K001
max. height = SHh + approx. 53 cm

Length: 55 cm
A9041401075-K002
max. height = SHh + approx. 63 cm
continuous height adjustment,
Safety push handle holder mounted on the back cross bar



Aluminum push handles, shifted back with side handle offset 145 mm

A9041401080-K* € 86.90

continuous height adjustment,
Safety push handle holder mounted on the back cross bar
Length: 55 cm
Up to SB 38 cm or body weight 60 kg



Aluminum central push handle Foldable handle area

A9041401086-K € 99.60

Height-adjustable and foldable in the handle area via toothed joint
Holder welded in the middle of the back crossbar
Up to SB 38 cm or body weight 50 kg
Not compatible with headrest holder
Attention: Not suitable for carrying the wheelchair - grip the backrest crossbar, to do it.



Central push handle for brake handles Foldable grip area

A9041401070-K € 137.80


closed, continuous height-adjustable, with foam rubber handle cover.
Up to SB 38 cm or body weight 50 kg
Not compatible with headrest holder
Only in conjunction with parking knee lever brake for accompanying person or drum brake



Aluminum push handle closed, can be folded down via toothed joint

A9041401069-K € 378.50

continuous height adjustment,
Push handle holder shifted back
Up to SB 38 cm or body weight 60 kg



Equipment without push handles

A9000401000



Backrest

Fixed backrest

A9041400020
with screwed aluminum side guards
Backrest angle 90° to the floor
Backrest cannot be folded



Folding backrest

A9041400010-K € 297.40
7-way angle adjustable,
with integrated plastic-coated aluminum side
parts



Back equipment

If not specified, the "Adaptable back with back padding" will be delivered.

Adaptable back with back padding

A9041802010-K
Cover: Black, breathable, open surface
Equipment with back cross bar



Ergonomic back shell aluminum from RH 24 cm

A9000802210-K € 428.70
with lumbar support
Colour: Anthracite matt
(If not specified under color selection)
If required, please select the back cushion
separately



Back cushion for ergonomic back shell

A9000802220-K € 211.90
Cushion thickness: 2 cm
Cover:
 Black, breathable, open surface
 Black, breathable, closed surface
If no specification is made, the cover material
"Black, breathable, open surface" will be
delivered.



ERGO Back backrest bar for ergonomic back shell from RH 24 cm and SB 30 cm

A9000802019-K € 235.50
with recessed gripping area
Colour: Anthracite matt
(If not specified under color selection)
only compatible with central push handle
art. no. A9041401082
If required, please select the back shell and
back cushion separately



Equipment without back upholstery and back padding or back system

A9000802011 € - 53.00
Equipment with back cross bar




Seat upholstery

If not specified, the "Body Contour seat upholstery" will be delivered.


Body Contour

A9000801025-K



Seat plate aluminum

A9000801210-K € 284.20
Mounted between the frame tubes,
on bolted welded seat plate holders
Colour: Anthracite matt
(If not specified under color selection)



Equipment without seat upholstery

A9000801050



Colour Back equipment and Seat plate

Information on current color designs and availability at www.proactiv-gmbh.de
If not specified, the seat and backrest system will be delivered in the standard colour anthracite matt.

Anthracite RAL 7016 gloss

- A9000100179-995
- Back shell
- ERGO Back backrest bar
- Seat plate aluminum

in frame colour

- A9000100179-999 **€ 89.60**
- Back shell
- ERGO Back backrest bar
- Seat plate aluminum

For frame coating in two-colour design

Please specify if the required option should be coated in the base colour or the design colour

- A9000100179-999 **€ 89.60**
- Back shell
- basic colour
- design colour
- ERGO Back backrest bar
- basic colour
- design colour
- Seat plate aluminum
- basic colour
- design colour

If a different color is required, please select this from the wheelchair colour palette:

in colour from one-colour design:

- A9000100179-998 **€ 145.00**

Please enter the last 3 digits of the required colour:

A9000100170-

- Back shell
- ERGO Back backrest bar
- Seat plate aluminum

in colour from two-colour design:



- A9000100179-997 **€ 234.80**

Please enter the last 3 digits of the required colour:

A9000100170-

- Back shell
- ERGO Back backrest bar
- Seat plate aluminum

Seat cushion



Seat cushion in rectangular contour Inlet: Standard

- A9000805180-K* **€ 111.90**

Seat cushion with radius contour Inlet: Standard

- A9000800300-K* **€ 155.40**

Seat cushion height:

- 3 cm
- 4 cm
- 5 cm
- 6 cm

Seat cushion in rectangular contour Inlet: Anti-Decubitus

- up to 30x30 cm A9000805204-K* **€ 290.00**
- >30x30 up to 50x50 cm A9000805200-K* **€ 422.00**

Seat cushion with radius contour Inlet: Anti-Decubitus

- up to 30x30 cm A9000800311-K* **€ 342.80**
- >30x30 up to 50x50 cm A9000800310-K* **€ 465.90**

Seat cushion height:

- 4 cm
- 6 cm

Cover:

- Black, breathable, open surface
- Black, breathable, closed surface

If no specification is made, the cover material "Black, breathable, open surface" will be delivered.

- Seat cushion is existing**

You can also find the seat cushions on our separate order sheet: "PRO ACTIV seat and backrest systems"

Clothing guard

Clothing guard plug-in

- A9041403110-K **€ 362.30**
removable



Aluminum clothing guard

Fixed to the side panel
Only in combination with fixed backrest

in black gloss

A9000403100-K € 236.10

in frame colour (basic colour with two-colour design)

A9000403101-K € 302.00

in design colour (only possible with frame colour in two-colour design)


A9000403102-K € 335.00

Armrests

Flip-down armrest with arm padding 25 cm

A9041403140-K € 297.60

Not in combination with clothing guard



Arm padding 25 cm length mounted on clothing guard

A9000403151 € 161.00

Only in combination with clothing guard
Depending on the clothing guard adjustment, the upper edge of the arm padding is mounted approx. 20 to 24 cm above the seat tube.



Axle plate

Standard equipment with 12 mm sockets/axles. Further axle plates on request.

2-slot axle plate

A9000201000*
Quick-release axle Ø 12 mm
Rear SH:
24" wheel: approx. 36 - 42 cm
22" wheel: approx. 33.5 - 39.5 cm
20" wheel: approx. 31 - 37 cm




1-slot axle plate

A9000201001*
Only available in black (no colour selection)
Only in combination with 10° camber
Rear SH can only be adjusted in some cases



Extended bike plate

A9000201030* € 306.00
+ 4° camber integrated
Adjustable track



Wheel camber

*Total width with tight-mounted handrim.

1° drive wheel camber

A9000202010*
Total width = approx. SB + 20 cm*

4° drive wheel camber

A9000202030*
Total width = approx. SB + 26 cm*

6° drive wheel camber

A9000202040*
Total width = approx. SB + 30 cm*

8° drive wheel camber

A9000202080* € 74.80
Total width = approx. SB + 33 cm*
Wheel plate mounted on the outside

10° drive wheel camber

A9000202050* € 132.90
Total width = approx. SB + 39 cm*
Wheel plate mounted on the outside

E-Adaptations

Preparation for the "e-motion" / „e-motion DuoDrive“ mounting

A9000201200 € 488.30
Equipment: e-motion torque support, 1/2" axle and socket length 50.8 mm (2"). PRO ACTIV and Alber recommend anti-tipping supports are installed in pairs. Please choose "Adaptations" under the selection or order Alber anti-tipping supports directly from Alber.

Preparation for the "e-fix" mounting

A9000201210 € 488.30
Equipment: e-fix torque support, 1/2" axle and socket length 50.8 mm (2"). PRO ACTIV and Alber recommend anti-tipping supports are installed in pairs. Please choose "Adaptations" under the selection or order Alber anti-tipping supports directly from Alber, as well as BUDDY mount.

Brake system

If not specified, the "Brake system push to lock with standard brake handle 20° outwards" will be delivered.

Brake system push to lock Brake handle standard 20° outwards

A9000700002



Brake system push to lock Long brake handle

A9000701010-002

Long brake handle 20° outwards

A9000701010-004



Brake lever with M8 adapter and plastic ball

A9000706001 € 162.40

only in combination with knee lever brake system, standard or long brake handle



Brake system push to lock Brake handle to flip-down

A9000706008 € 154.00

Length: 9 cm



Knee lever brake with one-hand operation

A9000707020-K € 380.90

via Bowden cable

Operation:

right side

left side



Parking knee lever brake for accompanying person

A9000701020 € 450.60

Push handle operation via Bowden cable.



Drum brake

A9000702011-K* € 460.20

Operating from the push handle
Total width + 3 cm (see wheel camber)
Please note:

- Only in combination with drive wheel aluminum
- 36-spoke rim wheel
- Not in combination with one-hand drive, e-fix or e-motion



Knee lever brake with lock

Standard brake handle A9000700214 € 154.00

Long brake handle A9000700215 € 154.00

Locking block for brake lever
Prevents braking during pushing by accompanying persons
Not possible in combination with knee lever brake with one-hand operation.



Footrest

If not specified, the "continuous footrest" will be delivered.

Continuous footrest

A9041601010-K*
with tubular resting bracket
height-adjustable
from USL 26 cm to USL 47 cm possible
If the footrest is wanted with a "rubber coated footplate", please select this under "Footrest accessories"



Footrest folding up to one side

A9041601300-K* € 168.80

from USL 25 cm:
with continuous tubular resting bracket
height adjustable +/- 3 cm
cannot be removed
from USL 29 cm:
can also be pivoted outwards
If the footrest is wanted with a "rubber coated footplate", please select this under "Footrest accessories"



Footrest folds up to the rear with spring locking mechanism

A9041602005-K* € 279.00


with continuous tubular resting bracket
from USL 32 cm
If the footrest is wanted with a "rubber coated footplate", please select this under "Footrest accessories"





<p>Divided footrest</p> <p><input type="checkbox"/> A9041603200-K € 327.00</p> <p>with rubber coated footplate supports height-adjustable can be pivoted outwards and folded up to the side cannot be removed from USL 29 cm and SB 26 cm possible</p> 	<p>Footrest, positioned on the upper frame</p> <p><input checked="" type="checkbox"/> A9041606010-K</p> <p>height-adjustable from USL 20 cm to USL 25 cm possible, Please select the rubber coated footrest under "Footrest accessories".</p> 	
-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--

Footrest accessories

<p>Rubber coated footrest with 14 cm depth</p> <p><input type="checkbox"/> A9000610000-K € 122.20</p> <p>bolted, angle-adjustable</p> 	<p>Rubber coated footrest with 18 cm depth</p> <p><input type="checkbox"/> A9000610005-K € 136.00</p> <p>bolted, angle-adjustable</p> 	<p>Footrest side stop</p> <p><input type="checkbox"/> A9000610100 € 58.00</p> <p>adjustable only in combination with rubber coated footplate</p> 
------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

<p>Deflector rollers mounted on the footplate</p> <p><input type="checkbox"/> A9000600006 € 114.00</p> <p>only with rubber coated footrest</p> 	
-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--

Adaptations

<p>Anti-tipping support</p> <p><input type="checkbox"/> right side A9041301160-K001 € 193.20</p> <p><input type="checkbox"/> left side A9041301160-K002 € 193.20</p> <p>with tension spring, can be folded up removable using Quick-pin</p> 	<p>Outdoor wheel attachment with steering fork</p> <p><input type="checkbox"/> A9041900100 € 821.10</p> <p>Equipped with straight-running lock can be removed from the wheelchair via quick release lever with swivelling pneumatic tyre 12.1/2 x 2.1/4" caster wheel can be removed from the support tube via quick-release axle for maximum ground pressure at the steering attachment of 30 kg</p> 	
-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--

You will find our PRO ACTIV transport restraint systems on our separate "Transport restraint systems" order sheet



3" caster wheel

3" safety light-up wheel (80 mm)

A9000501110 € 68.70
puncture-proof



4" caster wheel

4" caster wheel „Standard“ solid rubber black

A9000501011-004
(82A), puncture-proof



4" Skate wheel (100 mm)

A9000501120 € 68.70
puncture-proof



4" safety light-up wheel (100 mm)

A9000501220 € 92.00
puncture-proof



4" Frog Legs caster wheel (100 mm)

A9000501230 € 206.30
with replaceable rubber 6 mm, puncture-proof
Colour: black, silver rim
Incl. hair guard spacer (5.3 mm)



5" caster wheel

5" caster wheel „Standard“ PU black

A9000501036-005
(82A), puncture-proof



5" caster wheel „Standard“ PU grey

A9000501031-005
(82A), puncture-proof



5" safety light-up wheel (125 mm)

A9000501221 € 92.00
puncture-proof



5" aluminum rim 3-spoke design “good shock absorption”

Silver
A9000501210 € 174.60

Black
A9000501211 € 174.60

with replaceable rubber 6 mm medium hard,
puncture-proof
Incl. hair guard spacer (3 mm)



5" Frog Legs caster wheel (125 mm)

A9000501232 € 206.30
with replaceable rubber 6 mm, puncture-proof
Colour: black, silver rim
Incl. hair guard spacer (5.3 mm)



6" caster wheel

6" caster wheel „Standard“ PU grey

A9000501010-006
(82A), puncture-proof



6" caster wheel “pneumatic tyre“

A9000501050 € 55.50
SH front min. 42 cm



Caster fork

If not specified, the “standard caster fork with screw axle” will be delivered.

Caster fork Standard

A9000500100-K001
with screwed axle

A9000500100-K002 € 162.50
with quick-release axle



Caster fork Heavy Duty for high loads

Extra weight approx. 100 - 150 gr.

A9000510111-K001 € 52.40
with screwed axle

A9000510111-K002 € 212.50
with quick-release axle



Suspension fork

inkl. 5" aluminum rim “good shock absorption”

A9000500160-K002 € 598.90
with quick-release axle



Caster fork spare wheel „Standard“

A9000500300-K € 288.10
Caster fork Standard with quick-release axles



Caster fork spare wheel „Heavy Duty“

A9000500301-K € 328.10
Caster fork Heavy Duty with quick-release axles



The spare caster fork does not have a caster wheel.
In case a caster wheel is requested, please put a cross on a second selection in the category “Caster wheels”.

Track distance

Caster wheel track extension



A9000500400* € 57.40
10 mm per side

A9000500401* € 57.40
15 mm per side

A9000500410* € 57.40
20 mm per side

Drive wheels

If not specified, the "drive wheel 36 spokes aluminum rim" will be delivered.

Drive wheel 18-spoke aluminum rim

- 20" (ETRTO 451-19)
A9000203052
- 22" (ETRTO 489-19)
A9000203062
- 24" (ETRTO 540-19)
A9000203122

with Liteline hub, radially spoked

Not in combination with:

- Puncture-proof tyres
- One-hand drive
- Drum brake
- Transport restraint connection
- Handbike

up to 70 kg body weight:

Sturz \ Gewicht	1°	4°	6°	8°	10°
bis 50 Kg	✓	✓	✓	✓	✓
51 - 60 Kg	✓	✓	✓	✓	
61 - 70 Kg	✓	✓			



Drive wheel 36-spoke aluminum rim

- 20" (ETRTO 451-19)
A9000203050
- 22" (ETRTO 489-19)
A9000203060
- 24" (ETRTO 540-19)
A9000203000*

as standard with Liteline hub (left image),
in combination with one-hand drive with
aluminum front hub (right image),
radially spoked



Drive wheel double rim 36-spokes

- 24" (ETRTO 540)
A9000203100 **€ 153.20**
with aluminum front hub, radially spoked
Use with higher loads



Spinergy wheel "Spox Sport" with 24 PBO spokes; rim, hub and spokes in black

- 24" (ETRTO 540)
A9000203550 **€ 761.30**



Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim and hub in black, spokes in white

- 24" (ETRTO 540)
A9000203554 **€ 929.30**

Recommendation: up to 70 kg body weight,
not in combination with handbike/traction
device use.



Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim, hub and spokes in black


- 24" (ETRTO 540)
A9000203552 **€ 929.30**

Recommendation: up to 70 kg body weight,
not in combination with handbike/traction
device use.



Handrims

x = dimension handrim attachment



Handrim attachment
general:

- x = 11 mm
- x = 19 mm

Handrim attachment
Curve + Gekko handrim:

- x = 10 mm
- x = 18 mm

If not specified, the handrim will be mounted narrow (11 mm / 10 mm).

Stainless steel Ø 20 mm

- 20" handrim for 20" wheel
A9000205010-001
- 22" handrim for 22" wheel
A9000205010-002
- 24" handrim for 24" wheel
A9000205010-004



Aluminum Ø 19 mm silver anodised

- 24" handrim for 24" wheel
A9000205019-004



Aluminum Ø 16 mm powder coated in frame colour

(basic colour with two-colour design)

- = Gloss surface
- = Matt surface

Note: Matt coating is not recommended due to the low grip.

- 20" handrim for 20" wheel
A9000205021-001 **€ 66.30**
- 22" handrim for 22" wheel
A9000205021-002 **€ 66.30**



Aluminum Ø 19 mm powder coated in frame colour

(basic colour with two-colour design)

- = Gloss surface
- = Matt surface

Note: Matt coating is not recommended due to the low grip.

- 20" handrim for 20" wheel
A9000205020-001 **€ 66.30**
- 22" handrim for 22" wheel
A9000205020-002 **€ 66.30**
- 24" handrim for 24" wheel
9000205020-004 **€ 66.30**



Aluminum Ø 16 mm polished and hard-anodised







- 20" handrim for 20" wheel
A9000205031-001 **€ 160.10**
- 22" handrim for 22" wheel
A9000205031-002 **€ 160.10**



Aluminum Ø 19 mm polished and hard-anodised

- 20" handrim for 20" wheel
A9000205030-001 **€ 160.10**
- 22" handrim for 22" wheel
A9000205030-002 **€ 160.10**
- 24" handrim for 24" wheel
A9000205030-004 **€ 160.10**



<p>Titan ø 19 mm</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205040-004 € 325.80</p> 	<p>Silicone handrim cover</p> <p><input type="checkbox"/> for 20" handrim A9000205100-001 € 226.00</p> <p><input type="checkbox"/> for 22" handrim A9000205100-002 € 226.00</p> <p><input type="checkbox"/> for 24" handrim A9000205100-004 € 226.00</p> <p>Surface without structure, compatible with handrim Ø 19 mm or Ø 20 mm.</p> 	<p>Ultra-Grip „very grippy“</p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000205070-001 € 390.70</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205070-002 € 390.70</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205070-004 € 390.70</p> 
<p>Ergo-Grip</p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000205080-001 € 396.80</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205080-002 € 396.80</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205080-004 € 396.80</p> <p>Attachment fitting to the rim</p> 	<p>Maxgripp with cell and natural rubber cover very good grip, soft</p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000200788 € 396.80</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000200789 € 396.80</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000200787 € 396.80</p> 	<p>Curve black anodised, ergonomic handrim profile</p> <p><input type="checkbox"/> 24" handrim for 24" wheel H 24 mm A9000205130-004 € 389.10</p> <p><input type="checkbox"/> 24" handrim for 24" wheel H 31 mm, size L A9000205140-004 € 389.10</p> 
<p>Curve silver anodised ergonomic handrim profile</p> <p><input type="checkbox"/> 24" handrim for 24" wheel H 31 mm, size L A9000205160-004 € 394.30</p> 	<p>Curve powder-coated, ergonomic handrim profile</p> <p><input type="checkbox"/> = Frame colour (gloss surface)</p> <p><input type="checkbox"/> = Black (gloss surface)</p> <p><input type="checkbox"/> 24" handrim for 24" wheel H 24 mm A9000205135-004 € 454.00</p> <p><input type="checkbox"/> 24" handrim for 24" wheel H 31 mm, size L A9000205145-004 € 454.00</p>  <p>Profilhöhe H: 24mm oder 31mm</p>	<p>Curve Grip, with Tetra-Grip coating, ergonomic handrim profile</p> <p><input type="checkbox"/> 24" handrim for 24" wheel H 24 mm A9000205155-004 € 505.50</p> 

<p>Quadro with Tetra-Grip coating</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205170-004</p> <p style="text-align: right;">€ 517.50</p> 	<p>Gekko black anodised</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205175-004</p> <p style="text-align: right;">€ 501.20</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p> 	<p>Gekko silver anodised</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205185-004</p> <p style="text-align: right;">€ 513.20</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p> 
---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

<p>Gekko S black anodised</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205195-004</p> <p style="text-align: right;">€ 501.20</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 24 mm</p> 	<p>Equipment without handrims</p> <p><input checked="" type="checkbox"/> A9000205990</p>	
----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------	--

Spoke guard

Finger- and spoke guard "transparent"
Not in combination with Spinergy drive wheels.



<p><input type="checkbox"/> for 20" wheel A9000206300-001</p> <p style="text-align: right;">€ 201.10</p>	<p><input type="checkbox"/> for 22" wheel A9000206300-002</p> <p style="text-align: right;">€ 201.10</p>	<p><input type="checkbox"/> for 24" wheel A9000206300-004</p> <p style="text-align: right;">€ 201.10</p>
---------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------

Tyres

If not specified, the "Schwalbe RightRun in black" tyres will be delivered.

Schwalbe „RightRun“ High-pressure tyres

- 22x1" (489-25)
A9000204001-002
 - 24x1" (540-25)
A9000204001-004
- Colour: Black
7.5 bar



Schwalbe „RightRun“ High-pressure tyres

- 20x1" (451-25)
A9000204000-001
 - 22x1" (489-25)
A9000204000-002
 - 24x1" (540-25)
A9000204000-004
- Colour: Grey, black wall
7.5 bar



Schwalbe „DownTown“ Profiled tyres

- 24x1" (540-25)
A9000204010-004
- 7.5 bar



<p>Schwalbe "MarathonPlus Evolution" with "smart guard" puncture protection and reflector strips</p> <p><input type="checkbox"/> 22x1" (489-25) A9000204030-002 € 47.80</p> <p><input type="checkbox"/> 24x1" (540-25) A9000204031-004 € 47.80</p> <p>Colour: Black 10 bar</p> 	<p>Schwalbe „Marathon Plus Evolution“ with "smart guard" puncture protection</p> <p><input type="checkbox"/> 24x1" (540-25) A9000204030-014 € 47.80</p> <p>Colour: Black 10 bar</p> 	<p>Schwalbe „Marathon Plus“ with "smart guard" puncture protection</p> <p><input type="checkbox"/> 24x1" (540-25) A9000204030-004 € 47.80</p> <p>Colour: Grey, black wall 10 bar</p> 
<p>Schwalbe „One“</p> <p><input type="checkbox"/> 24x1" (540-25) A9000204123-004 € 90.30</p> <p>Colour: Black 9 bar</p> 	<p>Puncture-proof "solid rubber tyres"</p> <p><input type="checkbox"/> 22x1" (489-25) A9000204200-002 € 119.10</p> <p>Colour: Grey</p> <p><input type="checkbox"/> 24x1" (540-25) A9000204201-004 € 119.10</p> <p>Colour: Black Not in combination with the 18 spoke drive wheel, Spinergy and one-hand drive.</p> 	

Accessories

24" MTB spare wheel set

with double rim (ETRTO 507), aluminum front hub, quick-release axle and stainless steel handrim
Mountain bike tyres, selectively:

- A9000203410 € 682.90
24x1.75" Schwalbe "Marathon Plus"
Brake setting compatible with 24x1" wheel



- A9000203413 € 682.90
24x2.1" Schwalbe "SMART SAM Sport"
Brake setting not compatible with 24x1" wheel



NOTE: We charge a surcharge for the use of drum brakes

- A9000200004 € 221.10

Through axle with assistance, Large push button

- A9000200096 € 102.70
- Interlock dimension: 122 / 110 mm
Axle diameter: 12 mm
Not in conjunction with lockable drive wheel axle clip "Tetra Clip",
not in combination with one-hand drive




Lockable drive wheel axle clip "Tetra Clip"






- A9000206400 € 161.60
- Not in combination with Spinergy drive wheels,
one-hand drive, drum brake or Liteline hub



Hand air compressor with pressure gauge

- A9000901100 € 103.20
- up to 14 bar



<p>Length-adjustable safety belt with buckle</p> <p><input type="checkbox"/> A9000902100 € 78.30</p> 	<p>Adjustable safety belt with hook and fastener</p> <p><input type="checkbox"/> A9000902110 € 62.60</p> 	<p>Crutch holder with flexible nylon cup</p> <p><input type="checkbox"/> A9000902310* € 133.40</p> 
<p>SIGMA lighting set</p> <p><input type="checkbox"/> A9000903270 € 104.80</p> <p>Consisting of front light and rear lights, chargeable via Micro USB socket, incl. Micro USB charging cable Lighting set compliant with road traffic regulations</p> 	<p>Transport bag for drive wheels</p> <p><input type="checkbox"/> for 24" wheel A9000903530-004 € 75.50 Ø 76 cm x 2 cm</p> <p><input type="checkbox"/> for 26" wheel A9000903530-006 € 75.50 Ø 83 cm x 2 cm</p> 	<p>Theray table</p> <p><input type="checkbox"/> Width 350 mm (up to SB 37 cm) € 466.20 A9000902800-001</p> <p><input type="checkbox"/> Width 400 mm (SB 38 - 41 cm) € 466.20 A9000902800-002</p> <p><input type="checkbox"/> Width 420 mm (from SB 42 cm) € 466.20 A9000902800-003</p> <p>removable table top made of Lexan, transparent, with front edging, with 2 positioning stops Ø 25 mm (for adaptation to frame tube Ø 25 mm)</p> 
<p>Further accessories (such as handbikes, steering attachments, etc.) see separate order form or www.proactiv-gmbh.de.</p>		

Notes and remarks:

Location, Date

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

Signature End customer/dealer

- 18 - 2024-07-25 / The prices are plus statutory value-added tax and shipping costs

= Included in the basic price = Position selected