



BUDDY 4all Active easy running wheelchair

Maximum user weight limit 100 kg

A9047100000 € 3,830.00

Quotation date: _____
Consultation was carried out by
Staff at the medical supplies shop: _____

Order date: _____
according to GTC PRO ACTIV Reha-Technik GmbH

Medical suppliers / company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Shipping address: same as company address
 Other (see below)

Name: _____
Street: _____
Postcode / city: _____


Customer name / commission:

Surname, forename: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
D.O.B: _____
Paying organisation: _____

Invoicing address: same as company address
 other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Body size:



Body height: _____ cm

Weight: _____ cm

Lower leg length: _____ cm

(anatomic dimension: hollow of the knee to the foot sole)

Indication: _____

This information is necessary for checking. If the information is missing, a check of the relevant items cannot be made!

Basic size

Seat width SB	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	<input type="checkbox"/> 31
	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38					
	The following dimensions are optionally available subject to surcharge:				<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46

	Size M (Frame length 35 cm (to the bend) and frame height 21 cm) <i>not in conjunction with the straight front end</i>									
	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	<input type="checkbox"/> 31

Seat depth ST	Size L (Frame length 38 cm (to the bend) and frame height 21 cm)									
	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35
	Size XL (Frame length 43 cm (to the bend) and frame height 21 cm)									
	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40

Back height RH	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29
	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39
	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45				
	Optional with size XL:			<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48				

Lower leg length USL	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	
	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	
	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44									
	Optional with size L + XL:			<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47						

Seat height back SHh	for 20" drive wheel										
	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37				
	for 22" drive wheel										
	<input type="checkbox"/> 33,5	<input type="checkbox"/> 34,5	<input type="checkbox"/> 35,5	<input type="checkbox"/> 36,5	<input checked="" type="checkbox"/> 37,5	<input type="checkbox"/> 38,5	<input type="checkbox"/> 39,5				
for 24" drive wheel											
	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input checked="" type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42				

Seat height front SHv (for 5" caster wheel)	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47
---	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------

Tipping point KP (front side rear tube to axis middle drive wheel)	<input checked="" type="checkbox"/> 5	<input type="checkbox"/> 6	<input type="checkbox"/> 7	<input type="checkbox"/> 8	<input type="checkbox"/> 9	<input type="checkbox"/> 10
--	---------------------------------------	----------------------------	----------------------------	----------------------------	----------------------------	-----------------------------

Please cross on the appropriate dimension (all dimensions are specified in cm).

Underlined values are standard settings.

The values SHh and SHv may deviate a little from the specified values depending on the wheel size, wheel camber as well as the front and rear frame angle!

Frame colour - one-colour design

Information on current colour designs and example images at: www.proactiv-gmbh.de/farben.php

Frame colour: Selection from colour palette

Frame colour will be specified with order

A9000100170-000

Raw, blasted, not coated

A9000100170-020

Note: In combination with moisture e.g. hand perspiration, can affect dark aluminum discoloration.

Powder colours

White RAL 9016

gloss A9000100170-023

matt A9000100170-123 € 42.10

Brillant silver

gloss A9000100170-011

Anthracite RAL 7016

gloss A9000100170-001

matt A9000100170-101 € 42.10

Black RAL 9005

gloss A9000100170-003

matt A9000100170-019 € 42.10

Pure orange RAL 2004

gloss A9000100170-022

matt A9000100170-122 € 42.10

Pink magenta RAL 4010

gloss A9000100170-037

matt A9000100170-137 € 42.10

Pastel blue RAL 5024

gloss A9000100170-036

matt A9000100170-136 € 42.10

Lemon green

gloss A9000100170-005

matt A9000100170-105 € 42.10

Yellow RAL 1023

gloss A9000100170-046

matt A9000100170-146 € 42.10

Red RAL 3001

gloss A9000100170-047

matt A9000100170-147 € 42.10

Lagoon blue

gloss A9000100170-052

matt A9000100170-152 € 42.10

Athletic grey

gloss A9000100170-050

matt A9000100170-150 € 42.10

Muddy grey

gloss A9000100170-051

matt A9000100170-151 € 42.10

Powder colours „Dormant“

Deep black effect silver

gloss A9000100170-007

Sky blue

gloss A9000100170-017

Navy blue

gloss A9000100170-045

matt A9000100170-145 € 42.10

Electric green

gloss A9000100170-015

matt A9000100170-115 € 42.10

Violet

gloss A9000100170-044

Fire red

gloss A9000100170-018

matt A9000100170-118 € 42.10

Frame colour - One-colour design special colours “liquid coating”

€ 339.40

Storm (ice blue)

gloss A9000100170-500

matt A9000100170-501

Coffee

gloss A9000100170-520

matt A9000100170-521

Prairie (rose gold)

gloss A9000100170-504

matt A9000100170-505

Pearl (mother of pearl-white)

gloss A9000100170-524

matt A9000100170-525

Berry

gloss A9000100170-526

matt A9000100170-527

Blue ocean

gloss A9000100170-522

matt A9000100170-523

Emerald

gloss A9000100170-514

matt A9000100170-515

Mint

gloss A9000100170-528

matt A9000100170-529

Chili

gloss A9000100170-516

matt A9000100170-517

Frame colour - Two-colour design

Note: In the two-colour design, both colours are the same gloss level as the design colour.

Two-colour design with one base colour (on the rear frame) and one design colour (on the front frame) € 260.70

Please select a color from below:

Base colour (on the rear frame)

White RAL 9016

A9000100170-200

Brilliant silver

A9000100170-201

Anthracite RAL 7016

A9000100170-202

Black RAL 9005

A9000100170-203

Pastel blue RAL 5024

A9000100170-206

Lemon green

A9000100170-207

Yellow RAL 1023

A9000100170-208

Lagoon blue

A9000100170-212

Athletic grey

A9000100170-210

Muddy grey

A9000100170-211

Design colour (on the front frame)

Black RAL 9005

gloss A9000100170-259

matt A9000100170-260

Lemon green

gloss A9000100170-250

matt A9000100170-251

Chili

gloss A9000100170-253

matt A9000100170-254

Pure orange RAL 2004

gloss A9000100170-255

matt A9000100170-256

Coffee

gloss A9000100170-278

matt A9000100170-279

Mint

gloss A9000100170-286

matt A9000100170-287

Emerald

gloss A9000100170-274

matt A9000100170-275

Storm (ice blue)

gloss A9000100170-264

matt A9000100170-265

Prairie (rose gold)

gloss A9000100170-266

matt A9000100170-267

Pearl (mother of pearl-white)

gloss A9000100170-282

matt A9000100170-283

Berry

gloss A9000100170-284

matt A9000100170-285

Blue ocean

gloss A9000100170-280

matt A9000100170-281

Colour - Attachment parts

Note: If not specified, the caster fork and bearing block will be mounted in black.

Drive wheel: aluminum rim, hub and spokes

Black A9000100171-003

Caster fork and bearing block

Silver A9000100173-011

Black A9000100173-003

Colour design of attachment parts according to customer request:

Drive wheel: Aluminium rim and design inlet¹

Frame colour A9000100171-999 € 242.30

Select from one-colour design (different from the frame color) A9000100171-998 € 280.50

Colour: _____

A9000100170-

Caster wheel: 5" aluminium rim² (3-spoke design)

Frame colour A9000100174-999 € 73.00

Select from one-colour design (different from the frame color) A9000100174-998 € 84.60

Colour: _____

A9000100170-

Special colour requests (frame or attachment parts)

Frame colour for single purchasing

Selection from RAL colour palette

A9000100170-100 € 402.30

Gloss level

gloss

matt

colour name

RAL number

Note: attachment parts are **not** possible in design colors

¹ Only in conjunction with Liteline hub (with drive wheel aluminum rim with 18 or 36 spokes, single rim)


² Only in conjunction with caster wheel 5" aluminum rim 3-spoke design good shock absorption A900050121*, as well as with suspension fork 5-spoke design A9000500160*

Frames

Straight front end

(only with frame size L and XL)


A9047100010-K*
Knee angle 105°
Bearing block permanently welded



Front end abducted shape

3.5 cm / side

A9047100030-K*
Knee angle 105°
Bearing block permanently welded
Abduction starts at:
Size M: ST 26 cm
Size L: ST 29 cm
Size XL: ST 34 cm



Custom-made seat width dimension 39 - 46 cm

A9047100060 € 411.80
with inlet in the upper cross tubes,
lower cross tube double bolted,
side panel screw connection via the
reinforcement inlet

Push handles

Safety push handles back-positioned

Length: 45 cm
A9041401075-K001
max. height = SHh + approx. 53 cm

Length: 55 cm
A9041401075-K002
max. height = SHh + approx. 63 cm
continuous height adjustment,
Safety push handle holder mounted on the
back cross bar



Aluminum push handles, shifted back with side handle offset 145 mm

A9041401080-K* € 86.90
continuous height adjustment,
Safety push handle holder mounted on the
back cross bar
Length: 55 cm
Up to SB 38 cm or body weight 60 kg



Aluminum central push handle Foldable handle area

A9041401086-K € 99.60
Height-adjustable and foldable in the handle
area via toothed joint
Holder welded in the middle of the back
crossbar
Up to SB 38 cm or body weight 50 kg
Not compatible with headrest holder
Attention: Not suitable for carrying the
wheelchair - grip the backrest crossbar, to do it.



Central push handle for brake handles Foldable grip area

A9041401070-K € 137.80
closed, continuous height-adjustable,
with foam rubber handle cover.
Up to SB 38 cm or body weight 50 kg
Not compatible with headrest holder
Only in conjunction with parking knee lever
brake for accompanying person or drum brake




Aluminum push handle closed, can be folded down via toothed joint

A9041401069-K € 378.50
continuous height adjustment,
Push handle holder shifted back
Up to SB 38 cm or body weight 60 kg



Equipment without push handles

A9000401000



Backrest

Fixed backrest

A9047400020
with screwed aluminum side guards
Backrest angle 90° to the floor
Backrest cannot be folded



Folding backrest

A9000400023 € 297.40
7-way angle adjustable,
with integrated plastic-coated aluminum side
parts



Back equipment

If not specified, the "Adaptable back with back padding" will be delivered.

Adaptable back with back padding

A9041802010-K
Cover: Black, breathable, open surface
Equipment with back cross bar



Ergonomic back shell aluminum from RH 24 cm

A9000802210-K € 428.70
with lumbar support
Colour: Anthracite matt
(If not specified under color selection)
If required, please select the back cushion separately



Back cushion for ergonomic back shell

A9000802220-K € 211.90
Cushion thickness: 2 cm
Cover:
 Black, breathable, open surface
 Black, breathable, closed surface
If no specification is made, the cover material "Black, breathable, open surface" will be delivered.



ERGO Back backrest bar for ergonomic back shell from RH 24 cm and SB 30 cm

A9000802019-K € 235.50
with recessed gripping area
Colour: Anthracite matt
(If not specified under color selection)
only compatible with central push handle art. no. A9041401082
If required, please select the back shell and back cushion separately



Equipment without back upholstery and back padding or back system

A9000802011 € - 53.00
Equipment with back cross bar




Seat upholstery

If not specified, the "Body Contour seat upholstery" will be delivered.


Body Contour

A9000801025-K




Seat plate aluminum

A9000801210-K € 284.20
Mounted between the frame tubes, on bolted welded seat plate holders
Colour: Anthracite matt
(If not specified under color selection)



Equipment without seat upholstery

A9000801050



Colour Back equipment and Seat plate

Information on current color designs and availability at www.proactiv-gmbh.de

If not specified, the seat and backrest system will be delivered in the standard colour anthracite matt.

Anthracite RAL 7016 gloss



A9000100179-999
 Back shell
 ERGO Back backrest bar
 Seat plate aluminum

in frame colour

A9000100179-999 € 89.60
 Back shell
 ERGO Back backrest bar
 Seat plate aluminum

<p>For frame coating in two-colour design Please specify if the required option should be coated in the base colour or the design colour</p> <p><input type="checkbox"/> A9000100179-999 € 89.60</p> <p><input type="checkbox"/> <u>Back shell</u></p> <p><input type="checkbox"/> basic colour</p> <p><input type="checkbox"/> design colour</p> <p><input type="checkbox"/> <u>ERGO Back backrest bar</u></p> <p><input type="checkbox"/> basic colour</p> <p><input type="checkbox"/> design colour</p> <p><input type="checkbox"/> <u>Seat plate aluminum</u></p> <p><input type="checkbox"/> basic colour</p> <p><input type="checkbox"/> design colour</p>	<p>If a different color is required, please select this from the wheelchair colour palette:</p> <p>in colour from one-colour design:</p> <p><input type="checkbox"/> A9000100179-998 € 145.00</p> <p>Please enter the last 3 digits of the required colour:</p> <p>A9000100170- <input type="text"/> <input type="text"/> <input type="text"/></p> <p><input type="checkbox"/> Back shell</p> <p><input type="checkbox"/> ERGO Back backrest bar</p> <p><input type="checkbox"/> Seat plate aluminum</p>	<p>in colour from two-colour design:</p> <p><input type="checkbox"/> A9000100179-997 € 234.80</p> <p>Please enter the last 3 digits of the required colour:</p> <p>A9000100170- <input type="text"/> <input type="text"/> <input type="text"/></p> <p><input type="checkbox"/> Back shell</p> <p><input type="checkbox"/> ERGO Back backrest bar</p> <p><input type="checkbox"/> Seat plate aluminum</p>
--	---	--

Seat cushion

 <p>Seat cushion in rectangular contour Inlet: Standard</p> <p><input type="checkbox"/> A9000805180-K* € 111.90</p>	 <p>Seat cushion with radius contour Inlet: Standard</p> <p><input type="checkbox"/> A9000800300-K* € 155.40</p>	<p>Seat cushion height:</p> <p><input type="checkbox"/> 3 cm <input type="checkbox"/> 5 cm</p> <p><input type="checkbox"/> 4 cm <input type="checkbox"/> 6 cm</p>
<p>Seat cushion in rectangular contour Inlet: Anti-Decubitus</p> <p><input type="checkbox"/> up to 30x30 cm € 290.00 A9000805204-K*</p> <p><input type="checkbox"/> >30x30 up to 50x50 cm € 422.00 A9000805200-K*</p>	<p>Seat cushion with radius contour Inlet: Anti-Decubitus</p> <p><input type="checkbox"/> up to 30x30 cm € 342.80 A9000800311-K*</p> <p><input type="checkbox"/> >30x30 up to 50x50 cm € 465.90 A9000800310-K*</p>	<p>Seat cushion height:</p> <p><input type="checkbox"/> 4 cm</p> <p><input type="checkbox"/> 6 cm</p>

Cover:

Black, breathable, open surface

Black, breathable, closed surface

If no specification is made, the cover material "Black, breathable, open surface" will be delivered.


Seat cushion is existing

You can also find the seat cushions on our separate order sheet: "PRO ACTIV seat and backrest systems"

Clothing guard

<p>Aluminum clothing guard Fixed to the side panel Only in combination with fixed backrest</p> <p>in black gloss</p> <p><input type="checkbox"/> A9000403100-K € 236.10</p>	<p>in frame colour (basic colour with two-colour design)</p> <p><input type="checkbox"/> A9000403101-K € 302.00</p>	 <p>in design colour (only possible with frame colour in two-colour design)</p> <p><input type="checkbox"/> A9000403102-K € 335.00</p>
---	--	--

<p>Aluminum clothing guard removable via lock function Adjustable following wheel arch Only in combination with folding backrest</p> <p>in anthracite matt</p> <p><input type="checkbox"/> A9000403215-K € 236.10</p>	<p>in frame colour (basic colour with two-colour design)</p> <p><input type="checkbox"/> A9000403216-K € 302.00</p>	 <p>in design colour (only possible with frame colour in two-colour design)</p> <p><input type="checkbox"/> A9000403217-K € 335.00</p>
---	--	--

<p>Carbon clothing guard removable via lock function</p> <p><input type="checkbox"/> A9000400013-K € 547.30</p> <p>Adjustable following wheel arch Only in combination with folding backrest</p> 	
---	--

Knee lever brake system, brake lever integrated in the clothing guard see "Brake systems"

Armrests

Armrest, height adjustable, integrated in the clothing guard

Including aluminum clothing guard with arm padding 25 cm long
Removable via lock function
Only in combination with folding backrest



in anthracite matt

A9000400014-K € 553.00

in frame colour (basic colour with two-colour design)

A9000400015-K € 603.80

in design colour (only possible with frame colour in two-colour design)

A9000400016-K € 629.20

Armrest, height adjustable, integrated in the clothing guard

A9000400017-K € 836.90

Including carbon clothing guard with arm padding 25 cm long
Removable via lock function
Only in combination with folding backrest

Axle plate

Standard equipment with 12 mm sockets/axles. Further axle plates on request.

2-slot axle plate

A9000201000*
Quick-release axle Ø 12 mm
Rear SH:
24" wheel: approx. 36 - 42 cm
22" wheel: approx. 33.5 - 39.5 cm
20" wheel: approx. 31 - 37 cm




1-slot axle plate

A9000201001*
Only available in black
(no colour selection)
Only in combination with 10° camber
Rear SH can only be adjusted in some cases



Extended bike plate

A9000201030* € 306.00
+ 4° camber integrated
Adjustable track
Rear SH:
24" wheel: approx. 36 - 42 cm
22" wheel: approx. 33.5 - 39.5 cm
20" wheel: approx. 31 - 37 cm



Wheel camber

*Total width with tight-mounted handrim.

1° drive wheel camber

A9000202010*
Total width = approx. SB + 20 cm*

4° drive wheel camber

A9000202030*
Total width = approx. SB + 26 cm*

6° drive wheel camber

A9000202040*
Total width = approx. SB + 30 cm*

8° drive wheel camber

A9000202080* € 74.80
Total width = approx. SB + 33 cm*
Wheel plate mounted on the outside

10° drive wheel camber

A9000202050* € 132.90
Total width = approx. SB + 39 cm*
Wheel plate mounted on the outside

E-Adaptations

Preparation for the "e-motion" / „e-motion DuoDrive“ mounting

A9000201200 € 488.30

Equipment: e-motion torque support, 1/2" axle and socket length 50.8 mm (2"). PRO ACTIV and Alber recommend anti-tipping supports are installed in pairs. Please choose "Adaptations" under the selection or order Alber anti-tipping supports directly from Alber.

Preparation for the "e-fix" mounting

A9000201210 € 488.30

Equipment: e-fix torque support, 1/2" axle and socket length 50.8 mm (2"). PRO ACTIV and Alber recommend anti-tipping supports are installed in pairs. Please choose "Adaptations" under the selection or order Alber anti-tipping supports directly from Alber, as well as BUDDY mount.

Brake system

If not specified, the "Brake system push to lock with standard brake handle" will be delivered.

Recommendation: For KP 7 cm or lower, we recommend the „Brake system push to lock with brake handle 20° outwards“.

Brake system push to lock Standard brake handle

A9000701010-001



Brake system push to lock Long brake handle

A9000701010-002



Brake system push to lock Brake handle 20° outwards

Standard brake handle
A9000700002

Long brake handle
A9000701010-004



Brake lever with M8 adapter and plastic ball

A9000706001 € 162.40

only in combination with knee lever brake system, standard or long brake handle



Brake system push to lock Brake handle to flip-down

A9000706008 € 154.00

Length: 9 cm



Knee lever brake with one-hand operation

A9000707020-K € 380.90

via Bowden cable
Operation:

right side

left side



Knee lever brake system, brake lever integrated in the clothing guard

Brake lever mounted on the bottom frame tube

Operation via wire rope

Clothing guard cannot be removed.

Not possible with frame size M and tipping point >7 cm

in black gloss with fixed backrest
in anthracite matt with folding backrest

A9047701015-K € 624.30




in frame colour (basic colour with two-colour design)

A9047701016-K* € 689.60

in design colour (only possible with frame colour in two-colour design)

A9047701017-K* € 722.60



<p>Parking knee lever brake for accompanying person</p> <p><input type="checkbox"/> A9000701020 € 450.60</p> <p>Push handle operation via Bowden cable.</p> 	<p>Drum brake</p> <p><input type="checkbox"/> A9000702011-K* € 460.20</p> <p>Operating from the push handle Total width + 3 cm (see „wheel camber“) Please note: - Only in combination with drive wheel aluminum 36-spoke rim wheel - Not in combination with one-hand drive, e-fix or e-motion</p> 	<p>Knee lever brake with lock</p> <p><input type="checkbox"/> Standard brake handle A9000700214 € 154.00</p> <p><input type="checkbox"/> Long brake handle A9000700215 € 154.00</p> <p>Locking block for brake lever Prevents braking during pushing by accompanying persons Not possible in combination with knee lever brake with one-hand operation.</p> 
--	---	--

Footrest

If not specified, the "continuous footrest" will be delivered.

Continuous footrest

A9041601010-K*
with tubular resting bracket
height-adjustable
from USL 26 cm to USL 47 cm possible
If the footrest is wanted with a "rubber coated footplate", please select this under "Footrest accessories"



Footrest, positioned on the upper frame

A9041606010-K
height-adjustable
from USL 20 cm to USL 25 cm possible,
Please select the rubber coated footrest under "Footrest accessories".



Variable footrest positioning

A9047600001-K € 272.90
adjustable in height, depth and angle
incl. footplate support covered with anti-slip film



Footrest folds up to the rear with spring locking mechanism

A9047602005-K* € 337.10
with continuous tubular resting bracket
from USL 32 cm
including caster wheel track extension 20 mm per side
If the footrest is wanted with a "rubber coated footplate", please select this under "Footrest accessories"



Variable position foot support and folds up to the rear

A9047600002-K € 362.30
adjustable in height, depth and angle
incl. footplate support covered with anti-slip film



Footrest folding up to one side

A9041601300-K* € 168.80
from USL 25 cm:
with continuous tubular resting bracket
height adjustable +/- 3 cm
cannot be removed
from USL 29 cm:
can also be pivoted outwards
If the footrest is wanted with a "rubber coated footplate", please select this under "Footrest accessories"



Divided footrest

A9041603200-K € 327.00
with rubber coated footplate supports
height-adjustable
can be pivoted outwards and folded up to the side
cannot be removed
from USL 29 cm and SB 26 cm possible



Footrest accessories

Rubber coated footrest with 14 cm depth

- A9000610000-K € 122.20
bolted, angle-adjustable



Rubber coated footrest with 18 cm depth

- A9000610005-K € 136.00
bolted, angle-adjustable



Footrest side stop

- A9000610100 € 58.00
adjustable
only in combination with rubber coated footplate



Deflector rollers mounted on the footplate

- A9000600006 € 114.00



Adaptations

Anti-tipping support

- right side A9041301160-K001 € 193.20
 left side A9041301160-K002 € 193.20
with tension spring
can be folded up
removable using Quick-pin



Outdoor wheel attachment with steering fork

- A9041900100 € 821.10
Equipped with straight-running lock
can be removed from the wheelchair via quick release lever
with swivelling pneumatic tyre 12.1/2 x 2.1/4" caster wheel
can be removed from the support tube via quick-release axle
for maximum ground pressure at the steering attachment of 30 kg



You will find our PRO ACTIV transport restraint systems on our separate "Transport restraint systems" order sheet



3" caster wheel

3" safety light-up wheel (80 mm)

- A9000501110 € 68.70
puncture-proof



4" caster wheel

4" caster wheel „Standard“ solid rubber black

- A9000501011-004 (82A), puncture-proof



4" Skate wheel (100 mm)

- A9000501120 € 68.70
puncture-proof



4" safety light-up wheel (100 mm)

- A9000501220 € 92.00
puncture-proof



4" Frog Legs caster wheel (100 mm)

A9000501230 € 206.30

with replaceable rubber 6 mm, puncture-proof
Colour: black, silver rim
Incl. hair guard spacer (5.3 mm)



5" caster wheel

5" caster wheel „Standard“ PU black

A9000501036-005
(82A), puncture-proof



5" caster wheel „Standard“ PU grey

A9000501031-005
(82A), puncture-proof



5" safety light-up wheel (125 mm)

A9000501221 € 92.00
puncture-proof



5" aluminum rim 3-spoke design “good shock absorption”

Silver
A9000501210 € 174.60

Black
A9000501211 € 174.60

with replaceable rubber 6 mm medium hard,
puncture-proof
Incl. hair guard spacer (3 mm)



5" Frog Legs caster wheel (125 mm)

A9000501232 € 206.30
with replaceable rubber 6 mm, puncture-proof
Colour: black, silver rim
Incl. hair guard spacer (5.3 mm)



6" caster wheel

6" caster wheel „Standard“ PU grey

A9000501010-006
(82A), puncture-proof



6" caster wheel “pneumatic tyre”

A9000501050 € 55.50
SH front min. 42 cm



Caster fork

If not specified, the “standard caster fork with screw axle” will be delivered.

Caster fork Standard

A9000500100-K001
with screwed axle

A9000500100-K002 € 162.50
with quick-release axle



Caster fork Heavy Duty for high loads

Extra weight approx. 100 - 150 gr.

A9000510111-K001 € 52.40
with screwed axle

A9000510111-K002 € 212.50
with quick-release axle



Suspension fork

inkl. 5" aluminum rim “good shock absorption”

A9000500160-K002 € 598.90
with quick-release axle



Caster fork spare wheel „Standard“

A9000500300-K € 288.10

Caster fork Standard with quick-release axles



Caster fork spare wheel „Heavy Duty“

A9000500301-K € 328.10

Caster fork Heavy Duty with quick-release axles




The spare caster fork does not have a caster wheel.

In case a caster wheel is requested, please put a cross on a second selection in the category "Caster wheels".

Track distance

Caster wheel track extension



A9000500400* € 57.40

10 mm per side

A9000500401* € 57.40

15 mm per side

A9000500410* € 57.40

20 mm per side

For footbar that folds up to the rear with spring locking mechanism A9047602005* including

Drive wheels

If not specified, the "drive wheel 36 spokes aluminum rim" will be delivered.

Drive wheel 18-spoke aluminum rim

20" (ETRTO 451-19)
A9000203052

22" (ETRTO 489-19)
A9000203062

24" (ETRTO 540-19)
A9000203122

with Liteline hub, radially spoked

Not in combination with:

- Puncture-proof tyres
- One-hand drive
- Drum brake
- Transport restraint connection
- Handbike

up to 70 kg body weight:

Gewicht	Sturz				
	1°	4°	6°	8°	10°
bis 50 Kg	✓	✓	✓	✓	✓
51 - 60 Kg	✓	✓	✓	✓	
61 - 70 Kg	✓	✓			



Drive wheel 36-spoke aluminum rim

20" (ETRTO 451-19)
A9000203050

22" (ETRTO 489-19)
A9000203060

24" (ETRTO 540-19)
A9000203000*

as standard with Liteline hub (left image),
in combination with one-hand drive with
aluminum front hub (right image),
radially spoked



Drive wheel double rim 36-spokes

24" (ETRTO 540)
A9000203100 € 153.20


with aluminum front hub, radially spoked
Use with higher loads



<p>Spinergy wheel "Spox Sport" with 24 PBO spokes; rim, hub and spokes in black</p> <p><input type="checkbox"/> 24" (ETRTO 540) A9000203550 € 761.30</p> 	<p>Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim and hub in black, spokes in white</p> <p><input type="checkbox"/> 24" (ETRTO 540) A9000203554 € 929.30</p> <p>Recommendation: up to 70 kg body weight, not in combination with handbike/traction device use.</p> 	<p>Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim, hub and spokes in black</p> <p><input type="checkbox"/> 24" (ETRTO 540) A9000203552 € 929.30</p> <p>Recommendation: up to 70 kg body weight, not in combination with handbike/traction device use.</p> 
--	--	---

Handrims







x = dimension handrim attachment
















Handrim attachment general:
 x = 11 mm
 x = 19 mm



Handrim attachment Curve + Gekko handrim:
 x = 10 mm
 x = 18 mm

If not specified, the handrim will be mounted narrow (11 mm / 10 mm).




<p>Stainless steel Ø 20 mm</p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000205010-001</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205010-002</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205010-004</p> 	<p>Aluminum Ø 19 mm silver anodised</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205019-004</p> 	<p>Aluminum Ø 16 mm powder coated in frame colour (basic colour with two-colour design)</p> <p><input type="checkbox"/> = Gloss surface <input type="checkbox"/> = Matt surface</p> <p>Note: Matt coating is not recommended due to the low grip.</p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000205021-001 € 66.30</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205021-002 € 66.30</p> 
<p>Aluminum Ø 19 mm powder coated in frame colour (basic colour with two-colour design)</p> <p><input type="checkbox"/> = Gloss surface <input type="checkbox"/> = Matt surface</p> <p>Note: Matt coating is not recommended due to the low grip.</p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000205020-001 € 66.30</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205020-002 € 66.30</p> <p><input type="checkbox"/> 24" handrim for 24" wheel 9000205020-004 € 66.30</p> 	<p>Aluminum Ø 16 mm polished and hard-anodised</p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000205031-001 € 160.10</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205031-002 € 160.10</p> 	<p>Aluminum Ø 19 mm polished and hard-anodised</p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000205030-001 € 160.10</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205030-002 € 160.10</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205030-004 € 160.10</p> 

<p>Titan ø 19 mm</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205040-004 € 325.80</p> 	<p>Silicone handrim cover</p> <p><input type="checkbox"/> for 20" handrim A9000205100-001 € 226.00</p> <p><input type="checkbox"/> for 22" handrim A9000205100-002 € 226.00</p> <p><input type="checkbox"/> for 24" handrim A9000205100-004 € 226.00</p> <p>Surface without structure, compatible with handrim Ø 19 mm or Ø 20 mm.</p> 	<p>Ultra-Grip „very grippy“</p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000205070-001 € 390.70</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205070-002 € 390.70</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205070-004 € 390.70</p> 
<p>Ergo-Grip</p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000205080-001 € 396.80</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205080-002 € 396.80</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205080-004 € 396.80</p> <p>Attachment fitting to the rim</p> 	<p>Maxgripp with cell and natural rubber cover very good grip, soft</p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000200788 € 396.80</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000200789 € 396.80</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000200787 € 396.80</p> 	<p>Curve black anodised, ergonomic handrim profile</p> <p><input type="checkbox"/> 24" handrim for 24" wheel H 24 mm A9000205130-004 € 389.10</p> <p><input type="checkbox"/> 24" handrim for 24" wheel H 31 mm, size L A9000205140-004 € 389.10</p> 
<p>Curve silver anodised ergonomic handrim profile</p> <p><input type="checkbox"/> 24" handrim for 24" wheel H 31 mm, size L A9000205160-004 € 394.30</p> 	<p>Curve powder-coated, ergonomic handrim profile</p> <p><input type="checkbox"/> = Frame colour (gloss surface)</p> <p><input type="checkbox"/> = Black (gloss surface)</p> <p><input type="checkbox"/> 24" handrim for 24" wheel H 24 mm A9000205135-004 € 454.00</p> <p><input type="checkbox"/> 24" handrim for 24" wheel H 31 mm, size L A9000205145-004 € 454.00</p>  <p>Profilhöhe H: 24mm oder 31mm</p>	<p>Curve Grip, with Tetra-Grip coating, ergonomic handrim profile</p> <p><input type="checkbox"/> 24" handrim for 24" wheel H 24 mm A9000205155-004 € 505.50</p> 






<p>Quadro with Tetra-Grip coating</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205170-004</p> <p style="text-align: right;">€ 517.50</p> 	<p>Gekko black anodised</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205175-004</p> <p style="text-align: right;">€ 501.20</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p>  	<p>Gekko silver anodised</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205185-004</p> <p style="text-align: right;">€ 513.20</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p> 
---	--	---

<p>Gekko S black anodised</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205195-004</p> <p style="text-align: right;">€ 501.20</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 24 mm</p>  	<p>Equipment without handrims</p> <p><input checked="" type="checkbox"/> A9000205990</p>	
---	---	--

Spoke guard			
<p>Finger- and spoke guard "transparent" Not in combination with Spinergy drive wheels.</p>			
<p><input type="checkbox"/> for 20" wheel A9000206300-001</p> <p style="text-align: right;">€ 201.10</p>	<p><input type="checkbox"/> for 22" wheel A9000206300-002</p> <p style="text-align: right;">€ 201.10</p>	<p><input type="checkbox"/> for 24" wheel A9000206300-004</p> <p style="text-align: right;">€ 201.10</p>	

Tyres		
<p>If not specified, the "Schwalbe RightRun in black" tyres will be delivered.</p>		
<p>Schwalbe „RightRun“ High-pressure tyres</p> <p><input type="checkbox"/> 22x1" (489-25) A9000204001-002</p> <p><input type="checkbox"/> 24x1" (540-25) A9000204001-004</p> <p>Colour: Black 7.5 bar</p> 	<p>Schwalbe „RightRun“ High-pressure tyres</p> <p><input type="checkbox"/> 20x1" (451-25) A9000204000-001</p> <p><input type="checkbox"/> 22x1" (489-25) A9000204000-002</p> <p><input type="checkbox"/> 24x1" (540-25) A9000204000-004</p> <p>Colour: Grey, black wall 7.5 bar</p> 	<p>Schwalbe „DownTown“ Profiled tyres</p> <p><input checked="" type="checkbox"/> 24x1" (540-25) A9000204010-004</p> <p>7.5 bar</p> 

<p>Schwalbe „MarathonPlus Evolution“ with “smart guard” puncture protection and reflector strips</p> <p><input type="checkbox"/> 22x1" (489-25) A9000204030-002 € 47.80</p> <p><input type="checkbox"/> 24x1" (540-25) A9000204031-004 € 47.80</p> <p>Colour: Black 10 bar</p> 	<p>Schwalbe „Marathon Plus Evolution“ with “smart guard” puncture protection</p> <p><input type="checkbox"/> 24x1" (540-25) A9000204030-014 € 47.80</p> <p>Colour: Black 10 bar</p> 	<p>Schwalbe „Marathon Plus“ with “smart guard” puncture protection</p> <p><input type="checkbox"/> 24x1" (540-25) A9000204030-004 € 47.80</p> <p>Colour: Grey, black wall 10 bar</p> 
<p>Schwalbe „One“</p> <p><input type="checkbox"/> 24x1" (540-25) A9000204123-004 € 90.30</p> <p>Colour: Black 9 bar</p> 	<p>Puncture-proof “solid rubber tyres”</p> <p><input type="checkbox"/> 22x1" (489-25) A9000204200-002 € 119.10</p> <p>Colour: Grey</p> <p><input type="checkbox"/> 24x1" (540-25) A9000204201-004 € 119.10</p> <p>Colour: Black Not in combination with the 18 spoke drive wheel, Spinergy and one-hand drive.</p> 	
<p>Accessories</p>		
<p>24" MTB spare wheel set</p> <p>with double rim (ETRTO 507), aluminum front hub, quick-release axle and stainless steel handrim Mountain bike tyres, selectively:</p>		
<p><input type="checkbox"/> A9000203410 € 682.90</p> <p>24x1.75" Schwalbe "Marathon Plus" Brake setting compatible with 24x1" wheel</p> 	<p><input type="checkbox"/> A9000203413 € 682.90</p> <p>24x2.1" Schwalbe "SMART SAM Sport" Brake setting not compatible with 24x1" wheel</p> 	<p>NOTE: We charge a surcharge for the use of drum brakes</p> <p><input type="checkbox"/> A9000200004 € 221.10</p>
<p>Through axle with assistance, Large push button</p> <p><input type="checkbox"/> A9000200096 € 102.70</p> <p>Interlock dimension: 122 / 110 mm Axle diameter: 12 mm Not in conjunction with lockable drive wheel axle clip "Tetra Clip", not in combination with one-hand drive</p> 	<p>Hand air compressor with pressure gauge</p> <p><input type="checkbox"/> A9000901100 € 103.20</p> <p>up to 14 bar</p> 	<p>Length-adjustable safety belt with buckle</p> <p><input type="checkbox"/> A9000902100 € 78.30</p> 

<p>Adjustable safety belt with hook and fastener</p> <p><input type="checkbox"/> A9000902110 € 62.60</p> 	<p>Crutch holder with flexible nylon cup</p> <p><input type="checkbox"/> A9000902310* € 133.40</p> 	<p>SIGMA lighting set</p> <p><input type="checkbox"/> A9000903270 € 104.80</p> <p>Consisting of front light and rear lights, chargeable via Micro USB socket, incl. Micro USB charging cable Lighting set compliant with road traffic regulations</p> 
<p>Transport bag for drive wheels</p> <p><input type="checkbox"/> for 24" wheel A9000903530-004 € 75.50 Ø 76 cm x 2 cm</p> 	<p>Therapy table</p> <p><input type="checkbox"/> Width 350 mm (up to SB 37 cm) € 466.20 A9000902800-001</p> <p><input type="checkbox"/> Width 400 mm (SB 38 - 41 cm) € 466.20 A9000902800-002</p> <p><input type="checkbox"/> Width 420 mm (from SB 42 cm) € 466.20 A9000902800-003</p> <p>removable table top made of Lexan, transparent, with front edging, with 2 positioning stops Ø 25 mm (for adaptation to frame tube Ø 25 mm)</p> 	
<p>Further accessories (such as handbikes, steering attachments, etc.) see separate order form or www.proactiv-gmbh.de.</p>		

Notes and remarks:

Location, Date

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

Signature End customer/dealer