



Quotation date: _____
Consultation was carried out by
Staff at the medical supplies shop: _____

Order date: _____
according to GTC PRO ACTIV Reha-Technik GmbH

One-hand drive with double hand rim

Approved up to 120 kg load capacity

Medical suppliers / company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Invoicing address: same as company address
 other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Shipping address: same as company address
 Other (see below)


Name: _____
Street: _____
Postcode / city: _____

Customer name / commission:

Surname, forename: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
D.O.B: _____
Paying organisation: _____
Indication: _____

Version, attachment for new wheelchair orders

Please note: In the case of a wheelchair with one-hand drive, the wheelchair attachment parts are only available in **black**.



Wheel camber: 1°, 4°, 6°, 8°, 10°, 12° for TRAVELER, SPEEDY A1 / A2, BUDDY classic / 4all
 Wheel camber: 2° for SPEEDY 4all / 4you / 4all Ergo / 4you Ergo, SPEEDY 4teen

| | | |
|---|--|---|
| TRAVELER incl. removable drive shaft | SPEEDY 4all / 4you / 4all Ergo / 4you Ergo / SPEEDY 4teen | SPEEDY A1 / A2 / BUDDY classic / BUDDY 4all |
| <input type="checkbox"/> A9030208010* € 1,461.10 | <input type="checkbox"/> A9016208010* € 1,287.20 | <input type="checkbox"/> A9010208010* € 1,287.20 |

Please select the required operating side:

Operation: Right Left

Version, mounting for retrofitting

Retrofit kit for PRO ACTIV wheelchair

incl. stainless steel handrim, high-pressure smooth-running or high-pressure profiled tyres.
 For further options for an extra charge, see wheelchair order sheet.

| | | |
|--|--|---|
| <input type="checkbox"/> on SPEEDY A1 / A2 | on TRAVELER incl. removable drive shaft | on SPEEDY 4all / 4you / 4all Ergo / 4you Ergo / 4teen One-hand axle tube powder-coated |
| <input type="checkbox"/> on BUDDY classic / 4all | <input type="checkbox"/> E8000208200 € 2,336.60 | <input type="checkbox"/> E8016200100-K € 139.40 |
| <input type="checkbox"/> E8000208100 € 2,039.90 | | Colour: _____ Seat width SB: _____ cm Product number: _____ |

Handrim A

Attachment: tight x = 11 mm wide x = 19 mm

Handrim Stainless steel

for 20" wheel, Ø 452 mm
 A9000205010-101
 for 22" wheel, Ø 488 mm
 A9000205010-102
 for 24" and 25" wheel Ø 539 mm
 A9000205010-104
 for 26" wheel, Ø 590 mm
 A9000205010-106

Handrim Aluminum powder coated in frame colour (basic colour with two-colour design)

= Gloss surface
 = Matt surface
 Note: Matt coating is not recommended due to the low grip.

for 20" wheel, Ø 452 mm
 A9000205020-101 € 33.20
 for 22" wheel, Ø 488 mm
 A9000205020-102 € 33.20
 for 24" and 25" wheel, Ø 538 mm
 A9000205020-104 € 33.20
 for 26" wheel, Ø 590 mm
 A9000205020-106 € 33.20

Handrim Aluminum polished and hard-anodised

for 20" wheel, Ø 452 mm
 A9000205030-101 € 80.00
 for 22" wheel, Ø 488 mm
 A9000205030-102 € 80.00
 for 24" and 25" wheel, Ø 538 mm
 A9000205030-104 € 80.00
 for 26" wheel, Ø 590 mm
 A9000205030-106 € 80.00

Handrim Ultra-Grip

for 20" wheel, Ø 452 mm
 A9000205070-101 € 195.40
 for 22" wheel, Ø 486 mm
 A9000205070-102 € 195.40
 for 24" and 25" wheel, Ø 548 mm
 A9000205070-104 € 195.40
 for 26" wheel, Ø 587 mm
 A9000205070-106 € 195.40

Handrim B with 3-spoke mounting

Please select handrim B at least 2" smaller than handrim A

Attachment: tight y = 28 mm wide y = 38 mm

Handrim Stainless steel

18", Ø 398 mm
 A9000208010-100
 20", Ø 452 mm
 A9000208010-101
 22", Ø 488 mm
 A9000208010-102
 24", Ø 539 mm
 A9000208010-104

Handrim Aluminum powder coated in frame colour (basic colour with two-colour design)

= Gloss surface
 = Matt surface
 Note: Matt coating is not recommended due to the low grip.

18", Ø 397 mm
 A9000208020-100 € 33.20
 20", Ø 452 mm
 A9000208020-101 € 33.20
 22", Ø 488 mm
 A9000208020-102 € 33.20
 24", Ø 538 mm
 A9000208020-104 € 33.20

Handrim Aluminum polished and hard-anodised

18", Ø 397 mm
 A9000208030-100 € 80.00
 20", Ø 452 mm
 A9000208030-101 € 80.00
 22", Ø 488 mm
 A9000208030-102 € 80.00
 24", Ø 538 mm
 A9000208030-104 € 80.00

Handrim Ultra-Grip

18", Ø 397 mm
 A9000208070-100 € 195.40
 20", Ø 452 mm
 A9000208070-101 € 195.40
 22", Ø 486 mm
 A9000208070-102 € 195.40
 24", Ø 548 mm
 A9000208070-104 € 195.40

Handrim C

Attachment: tight x = 11 mm wide x = 19 mm

Handrim Stainless steel

for 20" wheel, Ø 448 mm
 A9000205010-101
for 22" wheel, Ø 488 mm
 A9000205010-102
for 24" and 25" wheel, Ø 538 mm
 A9000205010-104
for 26" wheel, Ø 585 mm
 A9000205010-106

Handrim Aluminum powder coated in frame colour (basic colour with two-colour design)

= Gloss surface
 = Matt surface
Note: Matt coating is not recommended due to the low grip.

for 20" wheel, Ø 448 mm
 A9000205020-101 € 33.20
for 22" wheel, Ø 488 mm
 A9000205020-102 € 33.20
for 24" and 25" wheel, Ø 538 mm
 A9000205020-104 € 33.20
for 26" wheel, Ø 585 mm
 A9000205020-106 € 33.20

Handrim Aluminum polished and hard-anodised

for 20" wheel, Ø 448 mm
 A9000205030-101 € 80.00
for 22" wheel, Ø 488 mm
 A9000205030-102 € 80.00
for 24" and 25" wheel, Ø 538 mm
 A9000205030-104 € 80.00
for 26" wheel, Ø 585 mm
 A9000205030-106 € 80.00

Handrim Ultra-Grip

for 20" wheel, Ø 448 mm
 A9000205070-101 € 195.40
for 22" wheel, Ø 488 mm
 A9000205070-102 € 195.40
for 24" and 25" wheel, Ø 538 mm
 A9000205070-104 € 195.40
for 26" wheel, Ø 585 mm
 A9000205070-106 € 195.40

Equipment without handrim

A9000205990

Other

Silicone handrim cover Surface without structure

for handrim 20", Ø 448 mm
 A9000205100-001 € 112.90/piece
for handrim 22", Ø 488 mm
 A9000205100-002 € 112.90/piece
for handrim 24", Ø 538 mm
 A9000205100-004 € 112.90/piece
for handrim 26", Ø 585 mm
 A9000205100-006 € 112.90/piece

Please select which handrim you would like to attach the handrim cover to:

- handrim A
 handrim B
 handrim C

Knee lever brake with one-hand operation via bowden cable

A9000707020-K € 380.90
 Operation right side
 Operation left side

Knee lever brake for accompanying person

Operation from the push handle via Bowden cable

for TRAVELER

A9000701020 € 450.60

for SPEEDY A1 / A2

A9000701020 € 450.60

for SPEEDY 4teen

A9000701020 € 450.60

for BUDDY classic / 4all

A9000701020 € 450.60

Location, Date

Signature End customer/dealer

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.