



**UNIQUE LIKE YOU –  
FLEXIBLE LIKE YOUR LIFE**

The **SPEEDY 4you Ergo** active wheelchair impresses with its unique frame with fine geometry. 13 chassis parameters can be freely configured, including knee angle, frame retraction, wheelbase and lower leg width at the top and bottom. This makes the frame unique and precisely tailored to you. Thanks to the integrated tipping point and seat height adjustment in 1 cm increments, you can flexibly adapt the **SPEEDY 4you Ergo** to your needs at any time, even if these change over time.

**SPEEDY**  
**4you Ergo**




**#Individuality**

With 13 selectable frame parameters, every frame is unique. The frame stem can include a frame retraction.



**#Reliability**

Sturdy and durable frame for maximum mobility and compatibility with auxiliary drives, handbikes and traction devices.



**#Independence**

Easy folding of the backrest and tool-free removal of assemblies enable independent transport.



**#Flexibility**


Integrated tipping point and seat height adjustment in 1 cm steps for optimal adaptation to individual needs.



## PERFECTLY TAILORED TO YOU.

The **SPEEDY 4you Ergo** is manufactured exactly according to your requests. You can choose from a wide range of equipment, including footrests, brakes, push handles, back and seat systems, clothing guards, armrests, as well as drive and caster wheels. Your wheelchair, like you need it.

### UNIQUE, SELECTABLE OPTIONS FROM PRO ACTIV:



The 7-way angle-adjustable and foldable backrest offers maximum flexibility for the ideal sitting position.



The footrest, which can be folded towards the back with a spring locking mechanism, is easy to handle and flexible in everyday life.



Manoeuvrable one-hand drive with double handrim for users with hemiplegia or restricted mobility on one side of the body.



Elegant integral parking brake with comfortable and safe operation.

## TECHNICAL SPECIFICATIONS

Seat width:	30 - 52 cm (in 1 cm steps)
Seat depth:	36 - 54 cm (in 1 cm steps)
Back height:	20 - 48 cm (in 1 cm steps)
Lower leg length:	37 - 50 cm (in 1 cm steps)
Rear seat height:	38 - 47 cm (in 1 cm steps)
Front seat height:	42 - 54 cm (in 1 cm steps)
Tipping point:	3 - 11 cm (in 1 cm steps)
Total product weight:	Starting at 9.7 kg
Max. standard load weight:	Up to 120 kg
Max. load weight for Heavy Duty:	Up to 150 kg appropriate individual dimensions upon request

Data subject to design changes.

## FOLLOW US ON SOCIAL MEDIA



**PROACTIV**<sup>®</sup>

PRO ACTIV Reha-Technik GmbH

Im Hofstätt 11  
72359 Dotternhausen  
info@proactiv-gmbh.de

[www.proactiv-gmbh.com](http://www.proactiv-gmbh.com)

Rev. 09/2024

Subject to technical modifications, printing errors & colour deviations.