



# PROACTIV®



## 30cm



# LIFT activ

### WITH THE LIFT ACTIV YOU CAN REACH MORE.

The **LIFT activ** with seat height adjustment makes you flexible in everyday life as well as in different professions and workplaces. You can easily reach machines, furniture, kitchen equipment, printers and shelves. Conversion measures are often no longer necessary. At only 14.5 kg, the LIFT activ impresses with its lightweight, dynamic design and easy transport.



### #Independence

With the LIFT activ, you can master everyday life with even more self-determination. When shopping, you can reach goods from the upper shelves without outside help and are independent.



### #AtEyeLevel

The infinitely variable seat height adjustment raises the seat by up to 30 cm. In this position, the wheelchair user is at eye level with a person standing about 1.75 m tall.



### #Flexibility

Your LIFT activ adapts to you: The tipping point can be adjusted afterwards to meet your needs.



### #Independence

Thanks to the weight optimization, the foldable backrest and the tool-free removal of assemblies, you can easily transport your LIFT activ.




## YOUR KEY TO MORE INDEPENDENCE.

Your **LIFT activ** will be manufactured individually for you, to suit your needs. You have a large selection of equipment – from footrests and brakes to push handles, seat and back systems as well as drive and caster wheels. So your LIFT activ is exactly as you need it.

### UNIQUE, SELECTABLE OPTIONS FROM PRO ACTIV:



The 7-way angle-adjustable and foldable backrest offers maximum flexibility for the ideal sitting position.



The backrest, which can be folded down at the lumbar joint, allows an optimal pack size for transport – quickly and easily.



The removable and swivelling frame stem with split footrest offers plenty of space for your transfer.



The footrest, which can be folded towards the back, is easy to handle and flexible in everyday life.

## TECHNICAL SPECIFICATIONS

Seat width: 35 - 46 cm (in 1 cm steps)

Seat depth: 35 - 48 cm (in 1 cm steps)

Back height: 20 - 48 cm (in 1 cm steps)

Lower leg length: 30 - 48 cm (in 1 cm steps)

Rear seat height: 41 - 43 cm (in 1 cm steps)


Front seat height: 45 - 51 cm (in 1 cm steps)

Tipping point: 6 - 10 cm (in 1 cm steps)

Total product weight: Starting at 14.5 kg

Max. standard load weight: Up to 100 kg

## FOLLOW US ON SOCIAL MEDIA



**PROACTIV**<sup>®</sup>

PRO ACTIV Reha-Technik GmbH

Im Hofstätt 11  
72359 Dotternhausen  
info@proactiv-gmbh.de

[www.proactiv-gmbh.com](http://www.proactiv-gmbh.com)